

## SPECIFICATION

Wakogel(R) 60N, 150~425  $\mu$ m

for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (150~425 $\mu$ m)	under investigation
Particle size (coarser than 425 $\mu$ m)	under investigation
Precipitation volume	under investigation
Solubility in ethyl acetate	to pass test
Loss on drying at 180°C	under investigation
pH (50g/l aqueous slurry,25°C)	under investigation
Flowing rate	under investigation
Separation power	to pass test
Activity	to pass test