

Code No.23-0278

SPECIFICATION

WakosilR 100C18
for Column Chromatography

| REQUIREMENT | SPECIFICATION |
|---------------------------------------|------------------------------|
| Appearance | White ~ nearly white, powder |
| Particle shape | Spherical |
| Particle size (63~212 μm) | min.70% |
| Loss on drying at 85°C | max.5.0% |
| Practical test (1) | to pass test |
| Practical test (2) | to pass test |
| Total carbon (C) | 13.0~18.0% |
| Total nitrogen (N) | max.0.2% |