

Code No.24-0046

## SPECIFICATION

## Xylene Cyanol FF

Guaranteed Reagent

| REQUIREMENT  | SPECIFICATION   |
|--|---|
| * Appearance<br>Water-insoluble matter<br>Ethanol-insoluble matter<br>Absorbance (5mg/l)<br>Loss on drying at 105°C<br>Sensitivity | Dark green ~ dark purple,<br>crystalline powder<br>max.0.5%<br><br>max.1.0%<br><br>min.0.50<br><br>max.5.0%<br><br>to pass test |

JIS K 8272-1994

(\*:Additional test performed by Wako but not required by JIS)