

Code No.26-0013

SPECIFICATION

Zinc

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|--|
| * Appearance Copper (Cu) Cadmium (Cd) Lead (Pb) Arsenic (As) Iron (Fe) | Bluish white or silver gray, mossy max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.007%(mass/mass) max.0.1 ppm(mass/mass) max.0.005%(mass/mass) |

JIS K 8012:2006

*:Additional test performed by Wako but not required by JIS
This product is mossy.