

SPECIFICATION

Zinc Acetate Dihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [Zn(CH ₃ COO) ₂ ·2H ₂ O] Solubility in dil.acetic acid Chloride (Cl) Nitrate (NO ₃) Sulfate (SO ₄) Sodium (Na) Potassium (K) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe) Ammonium (NH ₄)	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.1ppm(mass/mass) max.5ppm(mass/mass) max.0.001%(mass/mass)

JIS K 8356 : 2020(Amendment 1)(*:Additional test performed by Wako but not required by JIS)