

Code No.26-0167

## SPECIFICATION

Zinc Stearate, 1.5  $\mu\text{m}$ 

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, powder
Particle diameter (d50)	max.1.5 $\mu\text{m}$
Loss on drying at 105°C	max.0.5%
Acidity [as CH <sub>3</sub> (CH <sub>2</sub> ) <sub>16</sub> COOH]	max.0.5%
Assay (as Zn) (volumetric analysis)	10.0~11.0%
Assay (capillary GC)	under investigation