

Code No.29-3185

SPECIFICATION

Presep(R)-C Silica Gel for Sample Pretreatment

| REQUIREMENT | SPECIFICATION |
|---------------------------|---------------|
| Appearance | to pass test |
| Loss on drying at 180°C | max.5.0% |
| Weight of column packing | 0.65~0.75g |
| Permitted pressure (1) | to pass test |
| Permitted pressure (2) | to pass test |
| Leakage of column packing | to pass test |
| Recovery | 95~105% |
| State of column packing | to pass test |