

Code No.35-4572

SPECIFICATION

Strong Basic Anion Exchange Resin No.1 (2% Cross-linking, 50-100 Mesh, Cl Type)

| REQUIREMENT                   | SPECIFICATION |
|-------------------------------|---------------|
| Whole uncracked beads         | min.90%       |
| Total exchange capacity (Cl)  | min.0.7eq/L   |
| Water retention capacity (Cl) | 65～75%        |
| Particle size <180micron      | max.10%       |
| Particle size >500micron      | max.10%       |