

Code No.01-0104

SPECIFICATION

L-Alanine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (C ₃ H ₇ NO ₂)(after drying)	min.99.0%(mass/mass)
Solubility in water	to pass test
Specific rotation [α] _{D20} (after drying)	+14.0~+15.0°
Loss on drying at 105°C	max.0.2%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Chloride (Cl)	max.0.02%(mass/mass)
Sulfate (SO ₄)	max.0.02%(mass/mass)
Heavy metals (as Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH ₄)	max.0.02%(mass/mass)
Foreign amino acids	to pass test

JIS K 9101:2014

(*:Additional test performed by Wako but not required by JIS)