

Code No.01-0227

SPECIFICATION

4-Aminoantipyrine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₁₁ H ₁₃ N ₃ O) Solubility in water Solubility in ethanol Melting point Residue after ignition (as sulfate) Suitability for phenol determination	Pale yellow ~ pale yellowish brown, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test to pass test 107~110°C max.0.05%(mass/mass) to pass test

JIS K 8048:2021(*:Additional test performed by Wako but not required by JIS)