

Code No.01-0330

SPECIFICATION

Ammonium Dihydrogenphosphate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (NH ₄ H ₂ PO ₄)	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	4.0~4.5
Chloride (Cl)	max.5ppm(mass/mass)
Nitrate (NO ₃)	max.0.001%(mass/mass)
Sulfate (SO ₄)	max.0.005%(mass/mass)
Sodium (Na)	max.0.005%(mass/mass)
Potassium (K)	max.0.005%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)

JIS K 9006:2023(*:Additional test performed by Wako but not required by JIS)