

Code No.01-0449

## SPECIFICATION

## Antimony(III) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION	
* Appearance	White ~ slightly yellow, crystals ~ crystalline powder min.98.0%	
Assay		
Solubility in dil. hydrochloric acid		to pass test
Sulfur compounds (as SO <sub>4</sub> )		max.0.01%
Sodium (Na)		max.0.005%
Copper (Cu)		max.5ppm
Lead (Pb)		max.0.003%
Arsenic (As)		max.0.002%
Iron (Fe)		max.0.002%

JIS K 8400-1994(\*:Additional test performed by Wako but not required by JIS)