

Code No.01-0673

# SPECIFICATION

## Ammonium Sulfate for Enzyme Refining

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	5.0~6.0
Copper (Cu)	max.0.2ppm
Magnesium (Mg)	max.0.1ppm
Calcium (Ca)	max.0.2ppm
Zinc (Zn)	max.0.1ppm
Cadmium (Cd)	max.0.1ppm
Lead (Pb)	max.0.3ppm
Arsenic (As)	max.0.1ppm
Manganese (Mn)	max.0.05ppm
Iron (Fe)	max.0.2ppm
Nickel (Ni)	max.0.2ppm
Assay	min.99.5%