

# SPECIFICATION

## Acetonitrile

for SHIMADZU Protein Peptide Sequencer

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (200nm)	max.0.05
Absorbance (220nm)	max.0.02
Absorbance (225~400nm)	max.0.01
Water	max.0.02%
Aldehydes (as CH <sub>3</sub> CHO)	max.0.5ppm
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.1ppm
Assay (capillary GC)	min.99.8%