

Code No.01-1396

## SPECIFICATION

### Agarose HE (High electroendosmosis) for Electrophoresis

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, powder
Solubility in water	to pass test
Gel point	to pass test
Electroendosmosis (-Mr)	0.20~0.25
Loss on drying at 105°C	max.10.0%
Sulfate (SO <sub>4</sub> )	max.0.4%
Jelly strength (1.5%)	min.800g/cm <sup>2</sup>