

Code No.01-1512

SPECIFICATION

Albumin, from Bovine Serum, Low Salt for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown or slightly greenish yellow, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying over silica gel in vacuum	max.5.0%
pH of a 10g/L solution at 25°C	4.5~6.5
Heavy metals (as Pb)	max.0.002%
Sodium (Na)	max.0.5%
Total nitrogen (N)	min.15.0%
Assay (electrophoresis)	min.98%