

Code No.01-1911

SPECIFICATION

Acetic Acid
for HPLC

| REQUIREMENT | SPECIFICATION |
|-----------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Solubility in water | to pass test |
| Absorbance (1→4, 250nm) | max.0.50 |
| Absorbance (1→4, 254nm) | max.0.10 |
| Absorbance (1→4, 260nm) | max.0.02 |
| Absorbance (1→4, 270nm) | max.0.005 |
| Absorbance (1→4, 300~400nm) | max.0.003 |
| Water | max.0.3% |
| Fluorescence test | to pass test |
| Assay (HPLC) | min.99.5% |