

Code No.01-2103

## SPECIFICATION

## 1mol/L Ammonium Formate Solution

for HPLC

| REQUIREMENT            | SPECIFICATION          |
|------------------------|------------------------|
| Appearance             | Colorless clear liquid |
| Absorbance (240nm)     | max.0.40               |
| Absorbance (254nm)     | max.0.03               |
| Absorbance (270~400nm) | max.0.01               |
| pH at 25°C             | 6.3~6.7                |
| Fluorescence test      | to pass test           |
| Concentration (20°C)   | 0.95~1.05mol/L         |