

Code No.01-2257

SPECIFICATION

(E)-Asarone

for the Japanese Pharmacopoeia General Tests

| REQUIREMENT | SPECIFICATION |
|--|---------------------|
| Description (JP method) | White powder |
| Identification (JP method) | to pass test |
| Purity (JP method) | - |
| Related substances (JP method) | to pass test |
| * Loss on drying over silica gel in vacuum | under investigation |
| * Assay (HPLC) | min.98.0% |

This product is suited for (E)-Asarone of general tests on The Japanese Pharmacopoeia.(*:Additional test performed by Wako)