

Code No.01-2412

SPECIFICATION

Asarinin

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Specific rotation [α]D ₂₀ (c=0.25,CHCl ₃)	under investigation
* Loss on drying over silica gel in vacuum	under investigation
* Assay (HPLC)	min.98.0%

*:Additional test performed by Wako