

## SPECIFICATION

Amino Acids Mixture Standard Solution, Type B  
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	measured value
Practical test	to pass test
4-Aminobutyric Acid	0.2250~0.2750 $\mu$ mol/mL
2-Aminoethanol	0.2250~0.2750 $\mu$ mol/mL
Ammonium Chloride	0.2250~0.2750 $\mu$ mol/mL
5-Hydroxy-DL-lysine	0.2250~0.2750 $\mu$ mol/mL
L-Ornithine	0.2250~0.2750 $\mu$ mol/mL
L-Lysine	0.2250~0.2750 $\mu$ mol/mL
1-Methyl-L-histidine	0.2250~0.2750 $\mu$ mol/mL
L-Histidine	0.2250~0.2750 $\mu$ mol/mL
3-Methyl-L-histidine	0.2250~0.2750 $\mu$ mol/mL
L-Anserine	0.2250~0.2750 $\mu$ mol/mL
L-Carnosine	0.2250~0.2750 $\mu$ mol/mL
L-Arginine	0.2250~0.2750 $\mu$ mol/mL