

Code No.02-0272

## SPECIFICATION

(2-Bromoethyl)benzene

Wako Special Grade

| REQUIREMENT           | SPECIFICATION                             |
|-----------------------|-------------------------------------------|
| Appearance            | Colorless ~ slightly yellow, clear liquid |
| Density (20°C)        | 1.365~1.375g/mL                           |
| Refractive index 20°C | 1.553~1.559                               |
| Water                 | max.0.1%                                  |
| Assay (capillary GC)  | min.95.0%                                 |