

## SPECIFICATION

### Bone Morphogenetic Protein 4 (truncated), Human, recombinant, Animal-derived-free for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ $\mu$ g
Electrophoresis test	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source: E. coli expressed human bone morphogenetic protein 4 (truncated) This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.