

Code No.02-1950

SPECIFICATION

Benzene-d₆, 99.5%(High purity) for Nuclear Magnetic Resonance

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.005%
Degree of deuteration	min.99.5%
Assay (capillary GC)	min.99.9%
Impurities (NMR) (1H)	to pass test
Impurities (NMR) (13C)	to pass test
Impurities (NMR) (31P)	to pass test
Impurities (NMR) (19F)	to pass test