

Code No.03-0174

SPECIFICATION

Diammonium Cerium(IV) Nitrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Yellowish red, crystals ~ crystalline powder min.95.0%
Assay	
Solubility in water	
Chloride (Cl)	
Sulfate (SO ₄)	
Sodium (Na)	
Potassium (K)	
Copper (Cu)	
Lead (Pb)	
Iron (Fe)	
	to pass test
	max.0.001%
	max.0.005%
	max.0.05%
	max.0.05%
	max.0.002%
	max.0.002%
	max.0.002%

JIS K 8556-1995

(*:Additional test performed by Wako but not required by JIS)