

Code No.03-0213

# SPECIFICATION

## Charcoal, Activated for Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Black, shot ~ small granules
Particle size (coarser than 300 $\mu\text{m}$ )	max.50%
Particle size (300~63 $\mu\text{m}$ )	min.40%
Particle size (finer than 63 $\mu\text{m}$ )	max.10%
Loss on drying at 105°C	max.15.0%
Substances darkened by ammonia solution	to pass test
Separation power	to pass test