

Code No.03-0300

## SPECIFICATION

## Cholesterol

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystals ~ powder
Solubility in ethanol	to pass test
Melting point	146~149°C
Specific rotation [ $\alpha$ ] <sub>D25</sub> (c=2, 1,4-dioxane) (after drying)	-34~-39°
Loss on drying at 105°C	max.0.2%
Acidity	to pass test
Assay (capillary GC)	min.90.0%