

## SPECIFICATION

## Copper(I) Chloride

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CuCl) Solubility in acid Sulfate (SO <sub>4</sub> ) Nitrogen compounds (as N) Sodium (Na) Potassium (K) Copper(II) (Cu <sup>2+</sup> ) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe)	White ~ slightly green, crystalline powder min.95.0%(mass/mass) to pass test max.0.2%(mass/mass) max.0.005%(mass/mass) max.0.05%(mass/mass) max.0.01%(mass/mass) max.0.5%(mass/mass) max.0.01%(mass/mass) max.0.02%(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass)

JIS K 8138:2013(Amendment 1)

(\*:Additional test performed by Wako but not required by JIS)