

Code No.03-0444

## SPECIFICATION

## Copper(II) Sulfate

Extra Pure

REQUIREMENT	SPECIFICATION
* Appearance	Slightly grayish white ~ grayish white, powder
Assay (CuSO <sub>4</sub> )	min.97.5%(mass/mass)
Solubility in dil.sulfuric acid	to pass test
Chloride (Cl)	max.0.003%(mass/mass)
Sodium (Na)	max.0.1%(mass/mass)
Potassium (K)	max.0.1%(mass/mass)
Lead (Pb)	max.0.05%(mass/mass)
Iron (Fe)	max.0.03%(mass/mass)

JIS K 8984:2018\*:Additional test performed by Wako but not required by JIS