

SPECIFICATION

Copper(II) Sulfate

Extra Pure

| REQUIREMENT | SPECIFICATION |
|--|---|
| * Appearance Assay (CuSO ₄) Solubility in dil.sulfuric acid Chloride (Cl) Sodium (Na) Potassium (K) Lead (Pb) Iron (Fe) | Slightly grayish white ~ grayish white, powder min.97.5%(mass/mass) to pass test max.0.003%(mass/mass) max.0.1%(mass/mass) max.0.1%(mass/mass) max.0.05%(mass/mass) max.0.03%(mass/mass) |

JIS K 8984:2018*:Additional test performed by Wako but not required by JIS