

## SPECIFICATION

## Creatinine

Wako Special Grade

| REQUIREMENT                 | SPECIFICATION                           |
|-----------------------------|---|
| Appearance                  | White ~ nearly white, crystals ~ powder |
| Solubility in water         | to pass test                            |
| Chloride (Cl)               | max.0.05%                               |
| Sulfate (SO <sub>4</sub> )  | max.0.005%                              |
| Heavy metals (as Pb)        | max.0.001%                              |
| Iron (Fe)                   | max.0.001%                              |
| Ammonium (NH <sub>4</sub> ) | max.0.02%                               |
| Assay                       | min.99.0%                               |