

Code No.03-1060

SPECIFICATION

Clofibrate

for Biochemistry

| REQUIREMENT | SPECIFICATION |
|--|---|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 1.139~1.146g/mL |
| Refractive index 20°C | 1.500~1.505 |
| Water | max.0.2% |
| Acidity (as C ₁₀ H ₁₁ ClO ₃) | max.0.2% |
| p-Chlorophenol (C ₆ H ₄ OH) | max.0.003% |
| Assay (capillary GC) | min.98.0% |