

# SPECIFICATION

## Chloroform 300

for Pesticide Residue-PCB Analysis

| REQUIREMENT                                    | SPECIFICATION                    |
|--|----------------------------------|
| Appearance                                     | Colorless clear liquid           |
| Density (20°C)                                 | 1.471~1.485g/mL                  |
| Water  | max.0.05%                        |
| Residue after evaporation                      | max.5ppm                         |
| Free chlorine                                  | to pass test (Cl:max.about 1ppm) |
| Suitability for pesticide residue-PCB analysis | to pass test                     |
| Assay (capillary GC)                           | min.99.0%                        |
| Ethanol (stabilizer)                           | 0.4~0.9%                         |