

SPECIFICATION

Chloroform 5,000 for Pesticide Residue-PCB Analysis

| REQUIREMENT | SPECIFICATION |
|--|----------------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 1.471~1.485g/mL |
| Water | max.0.05% |
| Residue after evaporation | max.5ppm |
| Free chlorine | to pass test (Cl:max.about 1ppm) |
| Suitability for pesticide residue-PCB analysis | to pass test |
| Assay (capillary GC) | min.99.0% |
| Ethanol (stabilizer) | 0.4~0.9% |