

Code No.03-1864

## SPECIFICATION

### Chloroform 5,000

for Pesticide Residue-PCB Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	1.471~1.485g/mL
Water	max.0.05%
Residue after evaporation	max.5ppm
Free chlorine	to pass test (Cl:max.about 1ppm)
Suitability for pesticide residue-PCB analysis	to pass test
Assay (capillary GC)	min.99.0%
Ethanol (stabilizer)	0.4~0.9%