

Code No.03-1938

## SPECIFICATION

Copper(II) Sulfate Solution ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ :10w/v%)  
for Waste Water Test

| REQUIREMENT  | SPECIFICATION     |
|--|-------------------|
| Appearance   | Blue clear liquid |
| Suitability for phenol determination                         | to pass test      |
| Assay (as $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ) (20°C) | 9.5~10.5w/v%      |