

Code No.03-2033

SPECIFICATION

(E)-Chlorogenic Acid

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White powder
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Specific rotation [ $\alpha$ ]D20(c=0.4,H2O)	under investigation
* Water	under investigation
* Assay (HPLC)	min.98.0%

\*:Additional test performed by Wako