

Code No.03-2318

SPECIFICATION

Cyclopentyl Methyl Ether, Super Dehydrated, with Stabilizer for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ nearly colorless, clear liquid
Density (20°C)	0.860~0.866g/ml
Refractive index 20°C	1.418~1.423
Water	max.0.001%
Residue after evaporation	max.0.02%
Acidity (as CH ₃ COOH)	max.0.001%
Peroxides (as H ₂ O ₂)	max.0.001%
Assay (capillary GC)	min.99.0%