

## SPECIFICATION

## 1-Chlorobutane

for SHIMADZU Protein Peptide Sequencer

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (230nm)	max.0.25
Absorbance (240nm)	max.0.05
Absorbance (254nm)	max.0.02
Absorbance (260~400nm)	max.0.01
Water	max.0.02%
Aldehydes (as CH <sub>3</sub> CHO)	max.1ppm
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.1ppm
Suitability for amino acids sequence analysis	to pass test
Assay (capillary GC)	min.99.0%