

Code No.03-2546

SPECIFICATION

CultureSureR Albumin, from Bovine Serum
for Cell Culture

REQUIREMENT	SPECIFICATION
Appearance	White ~ greenish yellow or yellowish brown, crystals ~ powder
Solubility in water	to pass test
Loss on drying over silica gel in vacuum	under investigation
pH of a 100g/L solution at 25°C	6.5~7.5
Bacterial endotoxins	max.3.0EU/mg
Globulin	to pass test
Mycoplasma test	to pass test
Free fatty acids	to pass test
Protease	to pass test
Assay (electrophoresis)	min.98%
Viable count test	max.1000CFU/g