

Code No.04-0377

## SPECIFICATION

1,4-Dioxane

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Density (20°C)	1.028~1.038g/mL
Freezing point	10.0~12.0°C
Water	max.0.5%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.05%
Assay (capillary GC)	min.98.0%
Stabilizer (BHT)	1~5ppm