

Code No.04-0424

SPECIFICATION

2,2'-Bipyridyl

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly brown, crystals ~ crystalline powder
Assay (C ₁₀ H ₈ N ₂)	min.99.5%(mass/mass)
Solubility in ethanol	to pass test
Melting point	69.5~72.0°C
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Suitability for iron determination	to pass test

JIS K 8486:2007(*:Additional test performed by Wako but not required by JIS)