

Code No.04-1790

## SPECIFICATION

## 2-Diethylaminoethyl Methacrylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.916~0.923g/mL
Refractive index 20°C	1.441~1.448
Water	max.0.5%
Acidity [as CH <sub>2</sub> :C(CH <sub>3</sub> )COOH]	max.0.5%
Assay (capillary GC)	min.95.0%