

Code No.04-2261

## SPECIFICATION

## Dextran 200,000

for Leukocyte Separation

REQUIREMENT	SPECIFICATION
Appearance	White powder
Limiting viscosity (after drying)	0.38~0.42
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D20(c=6,H2O) (after drying)	+193~+201°
Loss on drying at 105°C	max.5.0%
Residue after ignition (as sulfate)	max.0.1%
Reducing sugar (as glucose)	max.0.4%
Average molecular weight	180000~210000