

Code No.04-2845

SPECIFICATION

Dichloromethane for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (CH ₂ Cl ₂) (GC)	min.99.5%(mass/mass)
Density (20°C)	1.315~1.330g/mL
Water	max.0.01%(mass/mass)
Residue after evaporation	max.0.0002%(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test
* Suitability for pesticide residue analysis (GC-FPD)	to pass test
* Methanol (stabilizer)	0.2~0.5%(mass/mass)

JIS K 8117:2015(*:Additional test performed by Wako but not required by JIS)