

## SPECIFICATION

Dicloxacillin Sodium Salt Monohydrate Standard  
for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION	
Appearance	White ~ pale yellow, crystalline powder ~ powder to pass test	
Solubility in water		
Specific rotation [ $\alpha$ ] <sub>D20</sub> (c=1,H <sub>2</sub> O) (calculated on the dehydrous basis) Water		+128~+143°
Assay (HPLC)		3.0~4.5%
	min.98.0%	