

SPECIFICATION

Dichloromethane, Super Dehydrated
for Organic Synthesis

| REQUIREMENT | SPECIFICATION |
|--------------------------------------|---------------------------------------|
| Appearance | Colorless clear liquid |
| Water | max.0.001% |
| Residue after evaporation | max.0.001% |
| Acidity (as HCl) | max.0.001% |
| Free chlorine | to pass test (Cl:max.about 1ppm) |
| Substances darkened by sulfuric acid | to pass test |
| Iodine-consuming substances | to pass test (as HCHO:max.about 4ppm) |
| Assay (capillary GC) | min.99.5% |
| 2-Methyl-2-butene (stabilizer) | 0.0005~0.005% |