

Code No.04-3237

SPECIFICATION

Dicyclanil Standard  
for High Performance Liquid Chromatography

| REQUIREMENT                              | SPECIFICATION  |
|--|--|
| Appearance                               | White ~ pale yellow,<br>crystalline powder ~ powder<br>to pass test<br><br>under investigation<br><br>under investigation<br><br>min.95.0% |
| Solubility in N,N -dimethylformamide     |  |
| Melting point                            |  |
| Loss on drying over silica gel in vacuum |  |
| Assay (HPLC)                             |  |