

Code No.04-3257

## SPECIFICATION

### Diethylene Glycol Ethyl Methyl Ether

Wako Special Grade

| REQUIREMENT           | SPECIFICATION                              |
|-----------------------|--------------------------------------------|
| Appearance            | Colorless ~ nearly colorless, clear liquid |
| Density (20°C)        | 0.918~0.928g/mL                            |
| Refractive index 20°C | 1.408~1.412                                |
| Water                 | max.0.1%                                   |
| Assay (capillary GC)  | min.98.0%                                  |