

Code No.04-3319

SPECIFICATION

Diisopropyl Ether, Super Dehydrated
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.720~0.730g/ml
Refractive index 20°C	1.366~1.371
Water	max.0.001%
Residue after evaporation	max.0.02%
Acidity (as CH ₃ CH ₂ COOH)	max.0.01%
Peroxide (as H ₂ O ₂)	max.0.001%
Assay (capillary GC)	min.99.0%